



# Development of an aerosol-compatible cell culture exposure system and its application to quantify cellular uptake of particles at the air-liquid interface

Jessica R. Murray, PhD

US EPA: Office of Research and Development  
Center for Public Health and Environmental Assessment  
Public Health and Integrated Toxicology Division  
AIT Webinar: November 29, 2022






## Disclaimer

The research described in this presentation has been reviewed by the Center for Public Health and Environmental Assessment, U.S. Environmental Protection Agency. Approval does not signify that the contents necessarily reflect the views and policies of the Agency, nor does the mention of trade names of commercial products constitute endorsement or recommendation for use.





# Methodologically Challenging Chemicals Require Advanced Exposure Methods



Over 10% of the Toxic Substances Control Act (TSCA) inventory includes volatile organic compounds (VOCs) and insoluble compounds which are incompatible with high-throughput screening.





**To address this challenge, we need to accomplish the following:**

1. Develop **ALI exposure technology** to include VOCs and insoluble chemicals in screening efforts
2. Create **analytical dosimetry methods** to quantify deposition and cellular uptake
3. Identify **appropriate human lung cell models and endpoints** to protect human health



The Inhalation Toxicology Facilities Branch (ITFB) developed the EPA's **Cell Culture Exposure Systems (CCES)** which permits dynamic exposure of human lung cells to VOCs at **air-liquid interface (ALI)**.

- Medium-throughput: 6 doses + 4 technical replicates within standard 24-well cell culture plate
  - Allows Benchmark Dose (BMD) modeling to estimate *in vitro* Points of Departure (PODs) for portal of entry effects
- Real-time sampling allows accurate exposure conditions to be reported throughout 2 h exposure



- Heated enclosure is key to maintaining 37°C and >80% RH throughout 2 h exposure condition
- No changes in viability or TEER observed after 2 h exposure in CCES




# BMD Values Share Similarities to TLV Rank Order Following VOC Exposures at Air-liquid Interface (ALI)

## Benchmark Dose Modeling Approaches for Volatile Organic Chemicals Using a Novel Air-Liquid Interface *In Vitro* Exposure System

Adam M Speen ✉, Jessica R Murray, Quentin Todd Krantz, David Davies, Paul Evansky, Joshua A Harrill, Logan J Everett, Joseph L Bundy, Lisa A Dailey, Jazzlyn Hill ... Show more

*Toxicological Sciences*, Volume 188, Issue 1, July 2022, Pages 88–107, <https://doi.org/10.1093/toxsci/kfac040>



### Benchmark Dose Analysis:

- HTTr TempO-Seq analysis at sub-cytotoxic concentrations
- Comparative to representative *in vivo* LOEL/NOEL values
- Within a magnitude of ACGIH occupational exposure TLVs

Chemical Name	BEAS-2B Median BMD (ppm)	HPBE Median BMD (ppm)	Representative LOEL (ppm)	Representative NOEL (ppm)	TLV (ppm)
Acrolein	0.586	--	0.25	NR	0.1
1-Bromopropane	2.246	N/A	62.5	250	0.1
Formaldehyde	N/A	--	2	1	0.3
1,3-Butadiene	13.979	--	625	200	10
Carbon Tetrachloride	9.563	N/A	20	5	10
Acetaldehyde	N/A	--	400	150	25
Trichloroethylene	44.842	28.148	50	25	50
Dichloromethane	142.127	226.73	8400	4200	100

<b>Exposure Regimen</b>	2 h exposure at ALI in 24-well format, endpoints analyzed 4 h later
<b>Endpoints</b>	Viability (ATP), n=2; Cytotoxicity (LDH), n=4; TempO-Seq (n=2)
<b>Biological Replicates</b>	Conducted over three days, n=3

Office of Pesticide Programs (OPP) and Office of Pollution Prevention and Toxics (OPPT) nominated the following list for further evaluation:

- Didecyl dimethyl ammonium chloride: *antiseptic/disinfectant*
- Polyhexamethylene guanidine-phosphate: *disinfectant*
- O-phenylphenol: *biocide used as preservative*
- Metribuzin: *herbicide*
- Tetramethrin: *insecticide*
- Indoxacarb: *pesticide*
- Naled: *insecticide*
- Oxamyl: *pesticide*
- Azoxystrobin: *pesticide & fungicide*
- Zinc pyrithione: *fungistatic & bacteriostatic*


Prioritized given widespread use during COVID-19 pandemic



Science.org "Does disinfecting surfaces really prevent the spread of coronavirus?"

**Must be generated as aerosols:** utilized a Blaustein Atomizer Module (BLAM) paired with syringe pump to generate liquid aerosols at high particle concentrations with a narrow particle size distribution


# Transport Physics and Deposition Mechanisms Differ Between VOCs and Particles



$$m_p \frac{dv_p}{dt} = \underbrace{-f_{ric}(\vec{v}_p - \vec{v}_f)}_{\text{aerosols}} + \underbrace{m_p \vec{g}}_{\text{nanoparticles/VOCs}} + \underbrace{\vec{F}(t)}_{\text{nanoparticles/VOCs}}$$

## Langevin Equation for Transport

- Particle Acceleration ( $\geq 0.5 \mu\text{m}$ )
- Gravitation forces ( $> 0.5 \mu\text{m}$ )
- Diffusion/Brownian Motion ( $\ll 0.5 \mu\text{m}$ )



Klinger-Strobel et al, *Expert Opin Drug Deliv.* 2015


## Convection-Diffusion Equation:

$$\underbrace{\frac{\partial C_p}{\partial t}}_{\text{Mass change in the system}} = \underbrace{-\nabla \cdot (C_p \vec{U}_{air})}_{\text{Transport deposition by convection}} + \underbrace{D \nabla^2 C_p}_{\text{Transport deposition by diffusion}} + \underbrace{-\nabla \cdot (C_p \vec{U}_{particle})}_{\text{Transport deposition by external forces}}$$

Deposition = Impaction      Diffusion      Sedimentation



## 1. Aerosolize Fluorescent Tracer



## 2. CAD + CFPD Modeling




### Navier-Stokes Equation for Incompressible Flow

$$\frac{d}{dt} \int_{\Omega} \rho d\Omega + \int_S \rho(\vec{V} \cdot \vec{n}_S) dS = 0$$

$$\frac{d}{dt} \int_{\Omega} (\rho \vec{V}) d\Omega + \int_S (\rho \vec{V})(\vec{V} \cdot \vec{n}_S) dS + \int_S (\vec{\tau} \cdot \vec{n}) dS = \int_S (-p \vec{n}) dS$$

- Computer Aided Design (CAD) utilized to create replicas of exposure system
- Computational Fluid-Particle Dynamics (CFPD) Modeling applied: Eulerian-Lagrangian approach
- *Limitations:* System components must be modeled separately to minimize computational expense




# CAD Models of CCES Dilution Manifold for CFPD Simulation




## CFD Boundary Conditions & Assumptions

Software	ANSYS Fluent
Dilution Inlet Flows	Define velocity to match CCES Operational Parameters
Delivery Outlets	Define negative velocity to match CCES Operational Parameters
Main Exhaust	Define pressure, $P = 0$
Wall	Constant temperature (37°C), “no slip” boundary condition
Turbulence Model	Laminar, $Re < 20$
Particle Movement	Discrete Phase Method, assumes particles $\leq 10\%$ of total flow



## 1. Aerosolize Fluorescein, MMAD: 1.3 $\mu\text{m}$




Fluorescein Deposition: Original VOC CCES




## 1. Aerosolize Fluorescein, MMAD: 1.3 $\mu\text{m}$



Fluorescein Deposition: Original VOC CCES






## 2. CFD simulation of particle distribution centered at MMAD: 1.3 $\mu\text{m}$ , $\rho$ : 1.602 g/mL



CFD-DPM Particle Trace: Original VOC CCES





# CAD and CFD Streamline Prototype Testing





- Goal: Half-Log Serial Dilution
- Generation 1
- Generation 1 w/ New Flows
- Generation 2
- Generation 3
- Generation 4-v1
- Generation 4-v2


# CFD Predicted Performance of New Dilution Manifold



# CFD Predicted Performance of New Dilution Manifold






Fluorescein vs. Half-Log Target



# Asymmetrical Aerosol Flow Leads to Flow Splitter Failure


CFPD: Particle Trace





Dose 1

Dose 2


Dose 3









# Fluorescein Deposition Patterns on Filters Confirms CFD Predictions





# CFD-DPM Modeling vs. Empirical Testing for Full System





## CFD + Discrete Phase Method

- DPM → Wall Film vs. Trap boundary conditions tested to estimate deposition
  - DPM impingement is based on Weber number, which does not consider electrostatic forces
- Very time intensive: 12-24 h+ per simulation




# CFD-DPM Modeling vs. Empirical Testing for Full System




## CFD + Discrete Phase Method

- DPM → Wall Film vs. Trap boundary conditions tested to estimate deposition
  - DPM impingement is based on Weber number, which does not consider electrostatic forces
- Very time intensive: 12-24 h+ per simulation




## Fluorescein Deposition on 16HBE Cells









**Fluorescein Deposition on 16HBE Cells**



	Total (ng/cm <sup>2</sup> /h)	St. Dev. (%)	Dilution Ratio
<b>Nozzle 1</b>	155.8	5%	
<b>Nozzle 2</b>	47.2	9%	0.30
<b>Nozzle 3</b>	14.4	5%	0.30


# New System Delivers Both Aerosols and VOCs



Modular, patent-pending design can be configured for aerosol or VOC dilution and delivery.


- Within the same system, drastically different operational conditions are required to deliver aerosols vs. VOCs
  - Aerosol operation: 5 mL min<sup>-1</sup> per well
  - VOC operation: 12.5 mL min<sup>-1</sup> per well
- Further work is needed to adapt the aerosol generation system to produce 6 doses of particles for a diverse list of test agents

- Cell-free collection methods are desirable as a high-throughput, low-cost method to quantify cell deposition to:
  1. Test improvements to aerosol generation system
  2. Quantify performance of ACCES (or other ALI exposure devices) for a variety of aerosols when fluorescence-based detection methods are not an option
- Literature search yielded a wide range of reported cell-free deposition methods for the Vitrocell 24/48, a similar perpendicular flow system




# Cell-free Options to Estimate Deposition

Geometry Preserved?	Height of trumpet/nozzle	Insert?	Collection Method	Citation
NO	Not reported	No	18.5 mL DPBS in base	Majeed et al., <i>Toxicology Letters</i> , 2014
NO	Not reported	Yes, 8 mm stainless steel	8 mm Cambridge filters	Zhang et al., <i>Toxicology In Vitro</i> , 2022
NO	Not reported	Yes, 8 mm stainless steel	125 $\mu$ L cell culture media	Verstraelen et al., <i>ALTEX</i> , 2021
NO	Not reported	Yes, 24-well ThinCert	100 $\mu$ L DPBS in insert	Steiner et al., <i>Toxicology In Vitro</i> , 2018
NO	Not reported	Yes, 24-well Transwell	100 $\mu$ L PBS in insert	Giralt et al., <i>Toxicology Letters</i> , 2020





- Cell-free collection methods that change geometry of the ALI system will impact deposition
- To preserve geometry near the air-liquid interface, we paired basolateral collection methods with a membrane-free transwell







	Fluorescein (ng/cm <sup>2</sup> /h)	Rhodamine (ng/cm <sup>2</sup> /h)
16HBE Cells	552.6	250.9
Water	1012.6	28.0
Media (1% FBS)	452.7	133.6
TF1000 Filter	262.4	523.5

- Cell-free controls failed to provide reliable estimate of cell deposition
  - Cell culture media was the best option for estimating fluorescein deposition, but underestimated rhodamine deposition by 50%
- Both particles were generated under identical conditions, but compound-specific deposition patterns were observed
  - Fluorescein (-), MMAD: 1.7  $\mu\text{m}$
  - Rhodamine (+), MMAD: 1.3  $\mu\text{m}$







### ACCES Deposition: 24-Well Format



	Fluorescein (ng/cm²/h)	Rhodamine (ng/cm²/h)
16HBE Cells	552.6	250.9
Water	1012.6	28.0
Media (1% FBS)	452.7	133.6
TF1000 Filter	262.4	523.5




### MedTec CelTox Deposition: 6-Well Format



	Fluorescein (ng/cm²/h)	Rhodamine (ng/cm²/h)
16HBE Cells	63.9	64.8
Water	491.9	2.7
Media (1% FBS)	79.0	9.6
TF1000 Filter	21.2	55.7

- Across multiple ALI exposure devices, we **cannot utilize cell-free options** to reliably estimate deposition without validation
- Sophisticated **cell extraction and analytical detection methods** are required to quantify cell deposition for a given aerosol




Perkins et al, *Environ Toxicol Chem*, 2019



Input to AOP constructs cannot be characterized by exposure concentration:  
 → *deposition and cellular uptake are dependent on exposure system and cell system.*

We aimed to determine whether fluorescent tracers could distinguish between deposition and cellular uptake.

# Using Fluorescent Tracers to Distinguish Between Total Deposition and Cellular Uptake



## Experimental Approach:


### ALI Testing Conditions

- BLAM used to generate liquid particles:
  - Fluorescein MMAD: 1.7  $\mu\text{m}$
  - Rhodamine MMAD: 1.3  $\mu\text{m}$
- Krypton-85 ( $^{85}\text{Kr}$ ) used as charge neutralizer
- Samples analyzed immediately after ALI exposure (2 h duration)

### Submerged Testing Conditions



- Tested 100  $\mu\text{L}$ , 50  $\mu\text{L}$ , and 10  $\mu\text{L}$ 
  1. Same total dose ( $\text{ng}/\text{cm}^2$ )
  2. Same concentration ( $\text{ng}/\text{mL}$ )
- Test agents dissolved in HBSS for direct-dosing
- Samples analyzed after 2 h to match ALI exposure duration

# Characterization of ALI vs. Submerged Exposures



	ALI Exposure	Submerged, or Direct Liquid Application
<b>Pros</b>	<ul style="list-style-type: none"> <li>• Most physiologically relevant</li> <li>• Direct cell-toxicant interaction</li> <li>• Compatible with both VOCs and particles</li> </ul>	<ul style="list-style-type: none"> <li>• Easier and higher-throughput, no complex equipment required</li> </ul>
<b>Cons</b>	<ul style="list-style-type: none"> <li>• ALI exposure equipment is complex to operate and maintain</li> <li>• Difficult to quantify delivery to cell surface</li> </ul>	<ul style="list-style-type: none"> <li>• Incompatible with VOCs</li> <li>• Unknown cellular uptake</li> <li>• Volume is not standardized across direct liquid application studies</li> <li>• Liquid application disrupts ALI conditions               <ul style="list-style-type: none"> <li>○ Measurable changes in TEER and baseline transcriptomics</li> </ul> </li> </ul>


# Characterization of ALI vs. Submerged Exposures




NAS, Using 21st Century Science to Improve Risk-Related Evaluations, 2017

“To best inform evidence integration for risk assessment, *in vitro* studies should determine the relevant internal cellular target dose rates (amount per unit time) that result in the observed responses” – Phalen et al., *Journal of Aerosol Science*, 2021.


Fluorescein: ALI vs. Submerged (909 ng/cm<sup>2</sup>)




Exposure	Dose (µg/cm <sup>2</sup> )	Conc. (µg/mL)	%basolateral	%uptake	%yield
ALI	~0.909	N/A	96%	4%	
100 µL	0.909	3	34%	6%	84%
50 µL	0.909	6	56%	4%	83%
10 µL	0.909	30	75%	6%	91%




Fluorescein: ALI vs. Submerged (909 ng/cm<sup>2</sup>)



Fluorescein: Same Concentration (3 µg/mL)



Fluorescein: Volume vs Translocation



ALI vs. Submerged, Same Total Dose

Exposure	Dose (µg/cm <sup>2</sup> )	Conc. (µg/mL)	%basolateral	%uptake	%yield
ALI	~0.909	N/A	96%	4%	
100 µL	0.909	3	<b>34%</b>	6%	84%
50 µL	0.909	6	<b>56%</b>	4%	83%
10 µL	0.909	30	<b>75%</b>	6%	91%


Same Concentration

Exposure	Dose (µg/cm <sup>2</sup> )	Conc. (µg/mL)	%basolateral	%uptake	%yield
100 µL	0.909	3	<b>36%</b>	3%	87%
50 µL	0.455	3	<b>56%</b>	4%	92%
10 µL	0.091	3	<b>72%</b>	4%	93%




# Cellular Uptake of Rhodamine is Volume-Dependent


Rhodamine: ALI vs. Submerged (909 ng/cm<sup>2</sup>)



Rhodamine: Same Concentration (3 µg/mL)



Rhodamine: Volume vs Cellular Uptake



ALI vs. Submerged, Same Total Dose


Exposure	Dose (µg/cm <sup>2</sup> )	Conc. (µg/mL)	%basolateral	%uptake	%yield
ALI	0.120*	N/A	0%	100%	
100 µL	0.909	3	1%	<b>30%</b>	102%
50 µL	0.909	6	1%	<b>44%</b>	96%
10 µL	0.909	30	4%	<b>58%</b>	78%

Same Concentration

Exposure	Dose (µg/cm <sup>2</sup> )	Conc. (µg/mL)	%basolateral	%uptake	%yield
100 µL	0.909	3	12%	<b>45%</b>	84%
50 µL	0.455	3	17%	<b>64%</b>	90%
10 µL	0.091	3	17%	<b>83%</b>	59%


# Mucus Retention Leads to Lower Cellular Uptake






**Mucus Retention of Rhodamine:  
ALI vs. Submerged**



Exposure	Dose (µg/cm <sup>2</sup> )	Conc. (µg/mL)	%mucus	%uptake	%yield
ALI	~0.303	N/A	<b>77%</b>	23%	
100 µL	0.303	1	<b>29%</b>	55%	67%
50 µL	0.303	2	<b>53%</b>	29%	71%
10 µL*	0.303	10	--	36%	86%

\*It was not technically possible to recover 10 µL of the apical solution without disturbing and aspirating mucus.







Exposure Method

Exposure Method

Fluorescein  
Rhodamine

Fluorescein  
Rhodamine



- Direct liquid application is often proposed as a time- and cost-effective alternative to ALI exposures, but variable apical volumes used across liquid application studies will directly impact cellular uptake:
  - Internal dose ranged from 40 – 360 ng/cm<sup>2</sup> for a single exposure concentration in 16HBE cells



## Conclusions: Lessons Learned When Including Aerosols in ALI Screening Efforts

- Careful **validation and characterization are required for each test agent**
  - An ALI system optimized for VOC delivery may not be appropriate for aerosols without significant modifications
- CAD and CFD Modeling were time- and cost-effective approaches to redesign our exposure system and optimize operational parameters
- Fluorescent tracers can be recovered in cell lysate to quantify cell deposition and can also be applied to validate CFD models
- Cell-free controls are rarely appropriate to estimate cell deposition
- Exposure Concentration  $\neq$  Deposition  $\neq$  Cellular Uptake
  - *This is especially important for mucus-producing cell lines!*
- Submerged exposure conditions are not comparable to ALI exposures, and differences in cellular uptake must be considered when designing these studies
- Further work is needed to translate ALI deposition to Human Equivalent Concentrations (HEC) to support *in vitro* to *in vivo* extrapolation (IVIVE)



# Acknowledgements

## Thank you to our CPHEA and CCTE collaborators:

*Engineers, Toxicologists, Data Scientists, and Support Staff*

- Mark Higuchi, PhD
  - Jason Weinstein, MS
  - Todd Krantz
  - Paul Evansky
  - Michael Monsees
  - Elise Carlsten
  - Wyatt Zander
  - Dave Davies
  - George Hudson
  - Earl Puckett
  - Lisa Dailey
  - Josh Harrill, PhD
  - Logan Everett, PhD
  - Joseph Bundy, PhD
  - Rusty Thomas, PhD
  - Annie Jarabek, PhD
  - Shaun McCullough, PhD
- MedTec Biolab Inc.
- Jose Zavala, PhD
- UNC School of Medicine
- Adam Speen, PhD

Know anyone who is graduating soon? We are hiring a new postdoc!  
ORISE opportunity: EPA-ORD-CPHEA-PHITD-2022-07 at [zintellect.com](https://zintellect.com)